

BSI Standards Publication

Sizing parameters of surgical valve prostheses: Requirements regarding the application of ISO 5840-2



National foreword

This Published Document is the UK implementation of ISO/PAS 7020:2023.

The UK participation in its preparation was entrusted to Technical Committee CH/150/2, Cardiovascular implants.

A list of organizations represented on this committee can be obtained on request to its committee manager.

Contractual and legal considerations

This publication has been prepared in good faith, however no representation, warranty, assurance or undertaking (express or implied) is or will be made, and no responsibility or liability is or will be accepted by BSI in relation to the adequacy, accuracy, completeness or reasonableness of this publication. All and any such responsibility and liability is expressly disclaimed to the full extent permitted by the law.

This publication is provided as is, and is to be used at the recipient's own risk.

The recipient is advised to consider seeking professional guidance with respect to its use of this publication.

This publication is not intended to constitute a contract. Users are responsible for its correct application.

This publication is not to be regarded as a British Standard.

© The British Standards Institution 2023 Published by BSI Standards Limited 2023

ISBN 978 0 539 18190 6

ICS 11.040.40

Compliance with a Published Document cannot confer immunity from legal obligations.

This Published Document was published under the authority of the Standards Policy and Strategy Committee on 31 May 2023.

Amendments/corrigenda issued since publication

Date Text affected

PHRI ICIV

ISO/PAS

This is a preview of "PD ISO/PAS 7020:2023". Click here to purchase the full version from the ANSI store.

SPECIFICATION

First edition 2023-05

Sizing parameters of surgical valve prostheses: Requirements regarding the application of ISO 5840-2

Désignation des dimensions des prothèses valvulaires chirurgicales: Exigences relatives à l'application de l'ISO 5840-2



PD ISO/PAS 7020:2023 **ISO/PAS 7020:2023(E)**

This is a preview of "PD ISO/PAS 7020:2023". Click here to purchase the full version from the ANSI store.



COPYRIGHT PROTECTED DOCUMENT

© ISO 2023

All rights reserved. Unless otherwise specified, or required in the context of its implementation, no part of this publication may be reproduced or utilized otherwise in any form or by any means, electronic or mechanical, including photocopying, or posting on the internet or an intranet, without prior written permission. Permission can be requested from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office CP 401 • Ch. de Blandonnet 8 CH-1214 Vernier, Geneva Phone: +41 22 749 01 11 Email: copyright@iso.org Website: www.iso.org

Published in Switzerland

Foreword Introduction			v				
				1	Scop	pe	1
				2	Normative references		1
3	Teri	ms and definitions	1				
4	Abb	previated terms	4				
5	Information required for the outer container labelling						
	5.2	Intended valve to be replaced	5				
	5.3	Inflow orifice diameter	5				
		5.3.1 General 5.3.2 Mechanical valves 5.3.3 Stented bioprosthetic valves	5				
	5.4	Effective orifice diameter 5.4.1 General	7 7				
	5.5	5.4.2 Steady flow method Valve housing external diameter					
	5.6	External sewing ring diameter					
6	Lab	elling format	9				
Rihl	iogran	shv	10				

Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

The procedures used to develop this document and those intended for its further maintenance are described in the ISO/IEC Directives, Part 1. In particular, the different approval criteria needed for the different types of ISO document should be noted. This document was drafted in accordance with the editorial rules of the ISO/IEC Directives, Part 2 (see www.iso.org/directives).

ISO draws attention to the possibility that the implementation of this document may involve the use of (a) patent(s). ISO takes no position concerning the evidence, validity or applicability of any claimed patent rights in respect thereof. As of the date of publication of this document, ISO had not received notice of (a) patent(s) which may be required to implement this document. However, implementers are cautioned that this may not represent the latest information, which may be obtained from the patent database available at www.iso.org/patents. ISO shall not be held responsible for identifying any or all such patent rights.

Any trade name used in this document is information given for the convenience of users and does not constitute an endorsement.

For an explanation of the voluntary nature of standards, the meaning of ISO specific terms and expressions related to conformity assessment, as well as information about ISO's adherence to the World Trade Organization (WTO) principles in the Technical Barriers to Trade (TBT), see www.iso.org/iso/foreword.html.

This document was prepared by Technical Committee ISO/TC 150, *Implants for surgery*, Subcommittee SC 2, *Cardiovascular implants and extracorporeal systems*.

Any feedback or questions on this document should be directed to the user's national standards body. A complete listing of these bodies can be found at www.iso.org/members.html.

Introduction

0.1 General

In the past, inconsistencies have been reported with the labelling and instructions for use associated with sizing parameters and sizing procedures for surgical replacement heart valves, specifically, mechanical and stented bioprosthetic valves. These inconsistencies have led to confusion among some users about which size valve to implant in a particular patient and have also led to challenges in comparing results (published or otherwise) from one valve model to another. A solution to the problem can be achieved by providing more complete and accurate sizing information to the clinicians, which will ultimately benefit the patients.

ISO 5840-2 identifies a number of sizing parameters that are required in the labelling (including on the unit box, see ISO 5840-2:2021, 6.3.3, and instructions for use) to inform the selection of a surgical heart valve prosthesis to be implanted in a specific patient. However, no guidance is offered in ISO 5840-2 on how these parameters should be obtained.

0.2 Clinical rationale for additional sizing information

Successful valve replacement therapy requires that an adequate size surgical heart valve substitute is used, based on patient body size and the native valve annulus size. An understanding of valve sizing parameters and appropriate choice of size is critical to post-procedure success since a valve substitute that is too small for the patient can result in prosthesis-patient mismatch. For aortic valve replacements, severe mismatch has been reported in 5 % to 15 % of patients. [1] Severe prosthesis-patient mismatch leads to increased early, mid-term and late mortality, especially if the left ventricular ejection fraction (LVEF) is reduced. [2]-[5] In the mid-term, it causes a higher incidence of heart failure [5] and limits left ventricular mass regression. [5] In the long term, it can also contribute to accelerated structural valve degeneration (SVD). [6] Patients with severe prosthesis-patient mismatch can require replacement of the valve substitute with another having a larger effective orifice area (EOA). However, re-intervention has significant risk of mortality and morbidity.

The best approach to prosthesis-patient mismatch is prevention. This requires the surgeon to have clear and accurate information about the sizing parameters and EOA of each valve substitute.

A surgical heart valve substitute is described by a labelled size given by the manufacturer, which is assumed to be broadly consistent with the size of the patient native valve annulus for which the valve is intended. Literature reviews^[1] and studies of haemodynamic function commonly compare valve substitutes by labelled size, but there can be major differences between the patient native valve annulus diameter and the labelled size of the valve substitute.^{[8]-[11]} Intraoperative sizing is further complicated by the need for aortic supra-annular valves to fit within the aortic sinus. The disparity between labelled size and actual size means that echocardiographic or clinical comparisons based on labelled size can be misleading.

The issue of valve sizing is a complex problem and is being addressed in a stepwise fashion. The working group revising ISO 5840-2 proposed a first step toward greater transparency by requiring additional information be added to the unit box, namely, internal orifice diameter and effective orifice diameter. Although this information does not necessarily inform the surgeon on whether the valve would fit in the patient's annulus, it helps to estimate the internal orifice available for blood flow and thus indirectly the EOA. It is not feasible to use clinically measured EOA's since sizing information must be available before a surgical heart valve substitute is released for use in patients. Indeed, it can take a number of years to gather sufficient echocardiographic data to confirm the clinical EOAs. Furthermore, the use of echocardiographic data to help avoid prosthesis-patient mismatch has been criticized because of variability in the measurements obtained in vivo. [12] In vitro steady flow data have less variability and allow meaningful comparison of every design and size of surgical heart valve substitute under the same flow conditions. This information can be used by the surgeon to choose a specific valve substitute type and size based on more standardized parameters than labelled valve size. It is anticipated that further steps toward a standardised approach to sizing will be addressed in subsequent editions of ISO 5840-2.

PD ISO/PAS 7020:2023 **ISO/PAS 7020:2023(E)**

This is a preview of "PD ISO/PAS 7020:2023". Click here to purchase the full version from the ANSI store.

This document provides further specifications to explain these two parameters (i.e. internal orifice diameter and effective orifice diameter) and other sizing parameters. This document also guides the manufacturer in selecting reproducible methods to obtain these parameters and the degree of accuracy required.

Sizing parameters of surgical valve prostheses: Requirements regarding the application of ISO 5840-2

1 Scope

This document describes in vitro methods of measurement of the sizing parameters for surgical valves (referring to mechanical and stented bioprosthetic valves only here and hereafter). It represents a consensus reached among manufacturers, independent bioengineers and clinicians, and is underpinned by interlaboratory studies.

This document relates to surgical heart valve prostheses and is intended to be used in conjunction with ISO 5840-1:2021 and ISO 5840-2:2021. Where noted, the requirements of this document clarify certain requirements of ISO 5840-1 and/or ISO 5840-2. Specific methodologies are included for flexible leaflet (bioprosthetic) and rigid (mechanical) valves. Sutureless valves, stentless valves and valved conduits are not included.

2 Normative references

The following documents are referred to in the text in such a way that some or all of their content constitutes requirements of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

ISO 5840-1:2021, Cardiovascular implants — Cardiac valve prostheses — Part 1: General requirements

ISO 5840-2:2021, Cardiovascular implants — Cardiac valve prostheses — Part 2: Surgically implanted heart valve substitutes

3 Terms and definitions

For the purposes of this document, the following terms and definitions apply.

ISO and IEC maintain terminology databases for use in standardization at the following addresses:

- ISO Online browsing platform: available at https://www.iso.org/obp
- IEC Electropedia: available at https://www.electropedia.org/

3.1 effective orifice area EOA

Α

orifice area that has been derived from flow and pressure or velocity data

Note 1 to entry: For in vitro testing, EOA is defined as:

$$A_{\text{eo}} = \frac{q_{V_{\text{RMS}}}}{51.6\sqrt{\frac{\Delta p}{\rho}}}$$

where

 $A_{\rm eo}$ is the effective orifice area, expressed in cm²;