

This is a preview of "ISO 24173:2009". [Click here to purchase the full version from the ANSI store.](#)

First edition
2009-09-01

Microbeam analysis — Guidelines for orientation measurement using electron backscatter diffraction

Analyse par microfaisceaux — Lignes directrices pour la mesure d'orientation par diffraction d'électrons rétrodiffusés



Reference number
ISO 24173:2009(E)

© ISO 2009

This is a preview of "ISO 24173:2009". [Click here to purchase the full version from the ANSI store.](#)

PDF disclaimer

This PDF file may contain embedded typefaces. In accordance with Adobe's licensing policy, this file may be printed or viewed but shall not be edited unless the typefaces which are embedded are licensed to and installed on the computer performing the editing. In downloading this file, parties accept therein the responsibility of not infringing Adobe's licensing policy. The ISO Central Secretariat accepts no liability in this area.

Adobe is a trademark of Adobe Systems Incorporated.

Details of the software products used to create this PDF file can be found in the General Info relative to the file; the PDF-creation parameters were optimized for printing. Every care has been taken to ensure that the file is suitable for use by ISO member bodies. In the unlikely event that a problem relating to it is found, please inform the Central Secretariat at the address given below.



COPYRIGHT PROTECTED DOCUMENT

© ISO 2009

All rights reserved. Unless otherwise specified, no part of this publication may be reproduced or utilized in any form or by any means, electronic or mechanical, including photocopying and microfilm, without permission in writing from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office
Case postale 56 • CH-1211 Geneva 20
Tel. + 41 22 749 01 11
Fax + 41 22 749 09 47
E-mail copyright@iso.org
Web www.iso.org

Published in Switzerland

This is a preview of "ISO 24173:2009". [Click here to purchase the full version from the ANSI store.](#)

Contents

Page

Foreword	iv
Introduction	v
1 Scope	1
2 Normative references	1
3 Terms and definitions	1
4 Equipment for EBSD	7
5 Operating conditions	8
6 Calibrations required for indexing of EBSPs	13
7 Analytical procedure	16
8 Measurement uncertainty	17
9 Reporting the results	18
Annex A (informative) Principle of EBSD	19
Annex B (normative) Specimen preparation for EBSD	20
Annex C (informative) Brief introduction to crystallography and EBSP indexing, and other information useful for EBSD	26
Bibliography	42

Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

International Standards are drafted in accordance with the rules given in the ISO/IEC Directives, Part 2.

The main task of technical committees is to prepare International Standards. Draft International Standards adopted by the technical committees are circulated to the member bodies for voting. Publication as an International Standard requires approval by at least 75 % of the member bodies casting a vote.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. ISO shall not be held responsible for identifying any or all such patent rights.

ISO 24173 was prepared by Technical Committee ISO/TC 202, *Microbeam analysis*.

This is a preview of "ISO 24173:2009". [Click here to purchase the full version from the ANSI store.](#)

Introduction

Electron backscatter diffraction (EBSD) is a technique that is used with a scanning electron microscope (SEM), a combined SEM-FIB (focussed-ion beam) microscope or an electron probe microanalyser (EPMA) to measure and map local crystallography in crystalline specimens [1],[2].

Electron backscatter patterns (EBSPs) are formed when a stationary electron beam strikes the surface of a steeply inclined specimen, which is usually tilted at $\approx 70^\circ$ from normal to the electron beam. EBSPs are imaged via an EBSD detector, which comprises a scintillator (such as a phosphor screen or a YAG single crystal) and a low-light-level camera (normally a charge-coupled device, CCD). Patterns are occasionally imaged directly on photographic film.

By analysing the EBSPs, it is possible to measure the orientation of the crystal lattice and, in some cases, to identify also the phase of the small volume of crystal under the electron beam. EBSD is a surface diffraction effect where the signal arises from a depth of just a few tens of nanometres, so careful specimen preparation is essential for successful application of the technique [3].

In a conventional SEM with a tungsten filament, a spatial resolution of about $0,25 \mu\text{m}$ can be achieved; however, with a field-emission gun SEM (FEG-SEM), the resolution limit is 10 nm to 50 nm, although the value is strongly dependent on both the material being examined and on the instrument operating parameters. Orientation measurements in test specimens can be carried out with an accuracy of $\approx 0,5^\circ$.

By scanning the electron beam over a region of the specimen surface whilst simultaneously acquiring and analysing EBSPs, it is possible to produce maps that show the spatial variation of orientation, phase, EBSP quality and other related measures. These data can be used for quantitative microstructural analysis to measure, for example, the average grain size (and in some cases the size distribution), the crystallographic texture (distribution of orientations) or the amount of boundaries with special characteristics (e.g. twin boundaries). EBSD can provide three-dimensional microstructural characterization by its use in combination with an accurate serial sectioning technique, such as focussed-ion beam milling [4].

It is strongly recommended that EBSD users be well acquainted with both the principles of crystallography and the various methods for representing orientations (both of which are described in the existing literature in this field) in order to make best use of the EBSD technique and the data produced [5],[6].