First edition 2010-09-15

# Nanotechnologies — Endotoxin test on nanomaterial samples for *in vitro* systems — *Limulus* amebocyte lysate (LAL) test

Nanotechnologies — Essai de détection d'endotoxines sur des échantillons de nanomatériaux pour des systèmes in vitro — Essai au lysat d'amébocyte de Limule (LAL)



#### **PDF disclaimer**

This PDF file may contain embedded typefaces. In accordance with Adobe's licensing policy, this file may be printed or viewed but shall not be edited unless the typefaces which are embedded are licensed to and installed on the computer performing the editing. In downloading this file, parties accept therein the responsibility of not infringing Adobe's licensing policy. The ISO Central Secretariat accepts no liability in this area.

Adobe is a trademark of Adobe Systems Incorporated.

Details of the software products used to create this PDF file can be found in the General Info relative to the file; the PDF-creation parameters were optimized for printing. Every care has been taken to ensure that the file is suitable for use by ISO member bodies. In the unlikely event that a problem relating to it is found, please inform the Central Secretariat at the address given below.



## **COPYRIGHT PROTECTED DOCUMENT**

#### © ISO 2010

All rights reserved. Unless otherwise specified, no part of this publication may be reproduced or utilized in any form or by any means, electronic or mechanical, including photocopying and microfilm, without permission in writing from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office Case postale 56 • CH-1211 Geneva 20 Tel. + 41 22 749 01 11 Fax + 41 22 749 09 47 E-mail copyright@iso.org Web www.iso.org Published in Switzerland

## Contents

Foroword	.,
	v
Introduction	V
1 Scope	1
2 Terms and definitions	1
3 Abbreviated terms	2
<ul> <li>4 Pre-test considerations</li></ul>	3 3 3 3
<ul> <li>5 Test sample</li> <li>5.1 Aqueous dispersion</li> <li>5.2 Aqueous extract</li> </ul>	3 3 3
<ul> <li>6 Preparation of test sample</li> <li>6.1 Dispersion method</li> <li>6.2 Extraction method</li> <li>6.3 Concentration</li> <li>6.4 Storage of test sample</li> <li>6.5 Laboratory environment</li> </ul>	3 3 4 4 4 4
<ul> <li>7 Test methods</li></ul>	5 5 5 6 6
<ul> <li>8 Assessment of results</li> <li>8.1 General</li> <li>8.2 Guidance on application of test</li> </ul>	6 6 7
9 Test report	7
Annex A (informative) Examples of potential interferences to LAL test	8
Annex B (informative) Gel-clot method	9
Annex C (informative) Endpoint photometric method1	3
Annex D (informative) Kinetic method1	6
Bibliography1	9

## Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

International Standards are drafted in accordance with the rules given in the ISO/IEC Directives, Part 2.

The main task of technical committees is to prepare International Standards. Draft International Standards adopted by the technical committees are circulated to the member bodies for voting. Publication as an International Standard requires approval by at least 75 % of the member bodies casting a vote.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. ISO shall not be held responsible for identifying any or all such patent rights.

ISO 29701 was prepared by Technical Committee ISO/TC 229, Nanotechnologies.

## Introduction

Endotoxins (lipopolysaccharides LPS) are part of the outer membrane of the cell wall of Gram-negative bacteria such as *E. coli, Salmonella, Shigella, Pseudomonas, Neisseria, Haemophilus.* Endotoxins can cause a variety of systemic reactions in mammals, including humans, such as fever, disseminated intravascular coagulation, hypotension, shock and death: the responses are mediated by production of various kinds of cytokines, activation of the complement cascade, activation of the coagulation cascade, etc. Endotoxins are present in the ordinary environment. Since most test samples of nanomaterials intended for *in vitro* and *in vivo* test systems require various preparation procedures, endotoxins might contaminate the test nanomaterials if the samples are prepared without special care.

For the purpose of toxicity screening or biocompatibility testing of nanomaterials, or mechanism studies on the possible toxicity induced by nanomaterials, various cell-based *in vitro* test systems and *in vivo* animal models are being developed and employed. In *in vitro* test systems, macrophages and other relevant mammalian cells are frequently used as the test cells especially for nanomaterials because they are primarily the responsible surveillance cells in the body. However, these cells are highly reactive to endotoxins; therefore it is difficult to distinguish the response to endotoxins from that to nanomaterials. Consequently, contamination by endotoxins would confound the result of tests *in vitro*.

Contamination by endotoxins of test samples may be reduced if appropriate precautions are followed in preparation of the test sample. Therefore the preliminary detection of endotoxins is required to minimize the contamination by endotoxins or confirm the insignificant levels of endotoxins in the test sample. It is also important to quantify endotoxin levels for the adequate interpretation of data obtained by *in vitro* biological test systems.

Since endotoxins may contaminate medical devices and medicines for parenteral use, quantitative and semiquantitative assay methods to test for endotoxins both *in vivo* and *in vitro* have been developed and used for regulatory purposes as well as laboratory standard operational procedures for nanomaterials (see Reference [6]). The bacterial endotoxin test using *Limulus* amebocyte lysate (LAL) reagent has been developed as an *in vitro* assay method to test for the presence of endotoxin contamination as an alternative to the pyrogenicity test using rabbits, and methods are described in the pharmacopoeia of many countries.

This International Standard provides considerations for the application of the LAL test to nanomaterial samples intended for *in vitro* biological tests.